



Scanning Electron Microscopy (SEM)

William Vanderlinde

IARPA

College Park, MD

William.e.vanderlinde@ugov.gov

Learning Objectives

- **To describe basic principles of Scanning Electron Microscopy techniques for microelectronic device failure analysis**

Outline

- **Fundamental theory of SEM imaging and analysis**
- **Practical tips for getting a good SEM image**
- **Recent developments and new features for SEM of semiconductors**

SEM vs. optical microscopy

SEM

High resolution (few nm)
Large depth of field
X-ray elemental analysis


Requires vacuum
Best on conductive samples
Poor TV rate imaging
Low contrast on defects
Difficult to navigate
May damage devices

Optical


Resolution $\sim \lambda/2 \sim 250$ nm
Very shallow depth of field
Color, phase contrast
Bright field / dark field

No vacuum
conductive or insulating
Live image
High contrast from defects
Easy to navigate
Dielectrics are transparent

SEM principles



Sample SEM image



200 μm ———

Mag = 50 x


Scanning principles

The image perspective is as if you are looking down the column

The detector provides an apparent “source of illumination” to the image

The image should always be viewed with the detector at the top


Deceptive images



1 μm —————

Mag = 30,000 x

Deceptive images




1 μm —————

Mag = 30,000 x


Deceptive images

**Detector at top
(correct image)**



1 μm ————— Mag = 30,000 x


Detector at bottom



1 μm ————— Mag = 30,000 x


These are images of the same object!

In-lens secondary electron detector




Secondary electron detectors

In-chamber detector



2 μm — Mag = 5,000 x


In-lens detector




2 μm — Mag = 5,000 x

In-lens detector: Less topography, reduced charging (?), better signal from deep holes

Magnification




Sample




TV Monitor

$$\text{Magnification} = \frac{L}{l}$$

To increase magnification, decrease the raster area



Maximum useful magnification



Beam diameter d when translated to the monitor has diameter $d \cdot M$

Maximum useful magnification

Image in sharp focus

$$d * M \leq 0.2 \text{ mm}$$





Image not in sharp focus


$$d * M > 0.2 \text{ mm}$$



$$M_{\max} = 0.2 \text{ mm} / d$$

For $d = 5 \text{ nm}$, the maximum useful mag is 40,000 x

Polaroid Mag v.s. Screen Mag



Polaroid Photo
3.5" x 4.5"
(9 cm x 11.5 cm)



Video Monitor

Digital Images



Computer Monitor

16 November 2009



What is the true mag?
Always use a micron marker

ISTFA 2009

17

Brightness


$$\text{Brightness} = \frac{\text{current}}{\text{area} \times \text{solid angle}} = \frac{4 i}{\pi^2 d^2 \alpha^2}$$

i = beam current

d = beam diameter

α = convergence angle


π = Pi ~ 3.14



SEM Cathode Comparison

Source:	<u>Tungsten</u>	<u>LaB₆</u>	<u>Schottky Field Emission</u>
Vacuum: (torr)	10^{-5}	10^{-7}	10^{-8}
Brightness: (A/cm ² ·sr)	10^{+5}	10^{+6}	10^{+8}
Resolution:	10 nm	5 nm	1 nm
Lifetime (hours)	40-100	200-1000	>1000

SEM column



Effect of lenses and apertures


$$\text{Brightness} = \frac{4 i}{\Pi^2 d^2 \alpha^2}$$

Lenses decrease d but increase α


Apertures decrease α but also decrease i

We wish to have small d , small α , large i



Electron Beam-Sample Interaction Products



Electron energy spectrum



Secondary electrons: topography contrast



Secondary electrons: material contrast



Secondary electron coefficient ~ 0.1 for most materials

Exceptions: Carbon ~ 0.05

Gold ~ 0.2



**Secondary electron coefficient strongly depends on
surface roughness, sample cleanliness, tilt angle**

Edge brightness




2 μm ——— Mag = 10,000 x

Backscatter: atomic number contrast



Secondary vs. Backscatter imaging



2 μm — Mag = 6,000 x



Secondary electron image



2 μm — Mag = 6,000 x


Backscattered electron image

Electron back-scatter diffraction (EBSD)




EBSD Pattern

Orientation Imaging Microscopy (OIM)




FIB image




OIM image

Ratchev, Carbonell, Ho, Bender, De Wolf,
Verlinden, Proceedings ISTFA 2002, p. 61-66.

Pincushion distortion





**SEM Grid
(true appearance)**





**500 μm — Mag = 35 x
Low mag distorted image**

Electron beam charging





$E_2 \sim 0.4 \text{ keV to } 4.0 \text{ keV}$

Charging examples



20 μm — Mag = 500x
Beam energy = 10 keV
Negative (bright) charging

Charging examples




10 μm ——— Mag = 2000x

Beam energy = 1.0 keV

Positive (dark) charging

Beam energies that reduce charging


Material	E ₂ (keV)
Polyimide	0.4
Photo Resist	0.55 – 0.70
PVC	1.65
Teflon	1.82
Glass Passivation	2.0
GaAs	2.6
Quartz	3.0
Alumina	4.2



20 μm – Mag = 333x

Beam energy = 2.5 keV

Passive Voltage Contrast Example




From S. Bothra et al., “A New Failure Mechanism by Corrosion of Tungsten in a Tungsten Plug Process”, Proceedings of the IRPS, 1998, pp. 150-156.

Slide courtesy Ed Cole

What about beam damage?

- “Raster burn” is often just sample charging which goes away when sample is vented to atmosphere.
- Large current in a small area can cause carbon deposition, especially if vacuum is poor.
- Large current on delicate samples (polymers) can melt or cause electrostatic discharge.
- If the electron beam penetrates to the gate, CMOS transistors may see threshold voltage shifts. This can generally be annealed out at 150 C for ~30 minutes.

“Raster Burn”




200 nm —

Mag = 60,000x

Low-K materials


Easily damaged by high energy electron beam



TEM

TEM

SEM




H. Bender, R.A. Donaton, ISTFA 2001

L. Li-Lung, ISTFA 2005

Sample Prep – Cross Sections

- Mechanical polish or FIB
- Automated fracture tools
- Staining and delineation



Typical Junction stain:

10:3:1 Acetic:Nitric:HF with bright light


Can add a few drops of copper sulfate solution

1 μm — Mag = 10,000x

Typical oxide delineation: buffered oxide etch

Sample Prep – RIE Deprocessing

- Reactive ion etching
- Selective to metal or dielectric
- Can be directional or isotropic
- May produce “grass” artifacts




1 μm —

Mag = 6,550x

SEM of reactive ion etched circuit

Sample Prep – Wet Deprocessing

- Wet chemical etching
- Very selective to materials
- Usually isotropic and will undercut layers
- Can be difficult to control



10 μm ————— Mag = 4,000x

SEM of circuit after HF chemical strip

Getting a good image

- Mounting the sample
- Selecting beam voltage
- Selecting working distance
- Column alignment and de-gaussing
- Reduce or eliminate charging
- Focus and Stigmatize the image
- Adjust the contrast and brightness

Sample mounting

Wish list: fast and easy, mechanically rigid, conductive, doesn't damage sample

- **Spring clips & screw mounts – fast & easy, rigid, but may damage sample**
- **Carbon or silver paint – rigid, safe, but slow**
- **Metal tape, double sticky carbon dots – fast & easy, but not rigid**

Selecting beam voltage


Maximum brightness ~ beam voltage

- High voltage produces higher brightness and smaller spot size for higher resolution
- But lower beam voltage is more sensitive to surface detail and can reduce charging
- Use the lowest beam voltage you can


Electron beam energy vs. range & spot size

Beam energy (keV)	Spot size (nm)	Range in Al (μm)
1	2.4	0.028
3.5	1.5	0.22
5	1.3	0.41
10	1.1	1.32
20	1.0	4.19
30	1.0	8.24

Range calculated from the Kanaya-Okayama formula



LEO 1550 FE “Gemini” Column Specs




**d = 1.0 nm
@ 20 kV**

**d = 2.3 nm
@ 1 kV**


Resolution vs. beam voltage for LEO 1550 FE

Beam Voltage vs. image detail



2 μm ——— Mag = 10,000 x


10 kV beam voltage



2 μm ——— Mag = 10,000 x


20 kV beam voltage

Voltage vs. image detail



2 μm ——— Mag = 10,000 x


3.5 kV beam voltage



2 μm ——— Mag = 10,000 x


1 kV beam voltage

Sample Tilt




- Tilting toward the detector increases the secondary electron emission coefficient and allows more electrons to be collected

Composing an image



2 μm ——— Mag = 10,000 x

No tilt




2 μm ——— Mag = 10,000 x

30 degrees tilt
45 degrees rotation


Short working distance

- **Short working distance decreases the beam spot size and allows for sharper high magnification images**
- **For many SEMs, a working distance of 2 mm to 10 mm is optimum for high resolution**
- **Long working distance improves depth of field and reduces pincushion or barrel distortion**

Achieving short working distance



OK



Better

Reducing electron beam charging

- Good electrical contact between the sample and the holder
- Use in-lens detector or variable pressure
- Increase tilt angle
- Reduce beam voltage
- Adjust beam current, raster area, raster rate
- Sputter coat sample with ~ 2 to 10 nm of Au/Pd, Cr, Pt, Ir, or C

Sputter coating samples







Pro:

- Sputter coating completely eliminates charging
- It increases electron emission from the surface
- It decreases the range of secondary electrons, thus increasing the surface sensitivity of the image

Con:

- You are imaging the coating, not the sample
- You may see artifacts from the coating grains
- You must complete all sample preparation such as plasma etching or chemical staining prior to sputter coating
- It may be more difficult to use energy dispersive x-ray analysis if you sputter coated with metal

Coatings at 2kV; 300kX



Sputter coating materials


- Au – very coarse grain size
- Au/Pd – coarse grain size
- Cr – fine grain size (but oxidizes in air)
- Ir – fine grain size (but requires high RF power, not compatible with some sputter coaters)
- Pt – fine grain size

[Au, Au/Pd, Pt coatings can be removed with aqua regia]


- Carbon – not as good for imaging (poor electron emissivity and range) but is x-ray transparent and can be removed with an oxygen plasma.
- Carbon sputters very slowly in a Ga ion beam and is a good definition layer between FIB deposited metals

Astigmatism


- Spherical asymmetry of a lens
- Corrected by additional lens coils
- Same procedure for SEM and FIB



Stigmatism




2 μm ——— Mag = 10,000 x



Under focus




2 μm ——— Mag = 10,000 x




Over focus

**Objective lens
(Focus knob)**


Correct focus



2 μm ——— Mag = 10,000 x




Correct focus




2 μm ——— Mag = 10,000 x

**Correct focus
with stigmatism
corrections**

Focus/Stig movie



Picking a focus object




Avoid lines in only one direction

Dust particle –
OK if not charging up

Orthogonal lines –
also OK


10 μm — Mag = 990 x

Brightness and contrast




2 μm ——— Mag = 10,000 x

Line scan




Too much contrast




2 μm ——— Mag = 10,000 x

Line scan




Not enough contrast

Line scan



2 μm ——— Mag = 10,000 x

Intensity ↑



Position ———→

Stereo Imaging in the SEM



Red/blue anaglyph




Stereoscope for SEM images

Modern methods of displaying stereo images:

- difficult to display to large audiences**
- difficult to archive**
- vertical height measurements are tedious**

Vertical Height Measurements



$$Z = (X_1 - X_2) / 2 \sin (\theta/2) \quad \text{Eqn. 1}$$


Where:

X_1 = distance in microns along the x-axis from the eucentric point to the feature in the first photo


X_2 = distance in microns along the x-axis from the eucentric point to the feature in the second photo

θ = tilt change between photos


A new approach: Digital Elevation Models (DEM)



0° tilt




5° tilt




DEM

**Image recognition software finds features,
calculates elevation, forms digital model**

A new approach: Digital Elevation Models (DEM)



DEM with SEM contrast




DEM with color-height scale

- **DEM with SEM contrast looks realistic**
- **DEM with color-height scale emphasizes topography**

DEM requires about 1 minute to calculate
DEM can be manipulated in real time

Profile Analysis




Profile analysis provides artificial surface profile and rapid measurements in X and Z

Ultra-high Resolution SEM

Requirements for ultra-high resolution SEM:



- (1) An electron beam finely focused to a small spot at the sample surface.**
- (2) Sufficient electron beam current to produce good image contrast.**
- (3) An imaging signal which originates very close to the electron beam impact site.**



(1) and (2) generally require high electron beam voltage, which causes problems obtaining (3).

Gold-on-Carbon Resolution Sample

Gold Islands




100 nm ——— Mag = 500,000 x


A highly reflective pattern (gold islands) on a strongly absorbing substrate (carbon) allows very high resolution imaging at high beam voltage.

STEM-in-SEM

STEM-in-SEM



— 100 nm Mag = 100,000x




— 100 nm Mag = 216,000x

High resolution is achieved on thin samples in the SEM using a specially designed sample holder.


Forward scattered electron imaging

Uncoated poly-silicon



50 nm — Mag = 200,000 x

Uncoated photo-resist



0.5 μm — Mag = 40,000 x

High resolution is achieved on bulk samples in the SEM using a sample holder designed to collect 30 keV electrons forward scattered at a high incident angle. (W. Vanderlinde, ISTFA 2003)

Short Courses

- **Lehigh Electron Microscopy School:**
<http://www.lehigh.edu/~inmatsci/shortcourses/Microcourses.html>
- **Maryland Practical Aspects of Electron Microscopy Short Course:**
<http://www.life.umd.edu/pasem/scanning.htm>
- **Pittcon Electron Microscopy and Microanalysis Short Course**
<http://www.micromaterialsresearch.com/Short%20Course.html>

References

- **Scanning Electron Microscopy and X-Ray Microanalysis : A Text for Biologists, Materials Scientists, and Geologists, by Joseph I. Goldstein, et al., 3rd Edition (2003).**
- **Scanning Electron Microscopy, X-ray Microanalysis, and Analytical Electron Microscopy: A Laboratory Workbook, by Charles E. Lyman, et al. (1990).**
- **“A review of wet etch formulations for silicon semiconductor failure analysis”, by Thomas W. Lee, Microelectronics Failure Analysis Desk Reference, 4th Edition.**

On the web

- <http://www.microscopy-online.com/>

Buyer's guides, lists of news groups and list servers.

- <http://www.ou.edu/research/electron/www-vl/>

Virtual Library of Microscopy. Thousands of links.

- <http://micro.magnet.fsu.edu/>

On-line virtual SEM. Also check out the “silicon zoo”.

- <http://mse.iastate.edu/microscopy/home.html>

SEM tutorials and picture galleries.

- <http://www.mos.org/sln/SEM/>

Museum of Science – a very basic introduction to the SEM.

SEM Materials and Supplies

- Ernest F. Fullam, Inc.

<http://www.fullam.com/>

- Ted Pella, Inc.

<http://www.tedpella.com/>

- M.E. Taylor Engineering, Inc.

<http://www.semsupplies.com/>

